



## VENTILATOR PORTABLE RESPIRATOR CORONAVIRUS COVID-19

**€7,990.00**

Portable Respirator Coronavirus COVID-19  
Ventilator for use in ambulance, hospital or  
private residence.

---

**SKU:** YSAV201P

**Categories:** [Bioide](#), [Machines](#), [Ventilators](#)

### PRODUCT DESCRIPTION

For use in ambulance, hospital or private residence.

Instrument classification: Class II

Warranty: 1 year

After-sale Service: Online Technical Support

model name: portable ventilator

Rate: 2~99bpm

Scope of application: adult, child

Gas source: O<sub>2</sub>, N<sub>2</sub>O: 0.27-0.55MPa

Power: AC 220V±10% 50Hz±1Hz

tidal volume range: 50~1500 ml

SIMV rate: 2~20bpm

Flow meter: O<sub>2</sub>: 0.1~1.0L/min 1.1-10L/min

Ventilation mode: IPPV, SIPPV, SIMV, PEEP, MANUAL, SIGH

Flow rate of rapid oxygen supply: 25-75L/min

